

# Direct Headspace Analysis of VOCs in Water by FT-MRR

## (Fourier Transform - Molecular Rotational Resonance)

Brent J. Harris <sup>a,b</sup>, Justin L. Neill <sup>a</sup>, Matt T. Muckle <sup>a</sup>, Robin L. Pulliam <sup>a</sup>, Dave A. McDaniel <sup>a</sup>, Roger L. Reynolds <sup>a</sup>, Brooks H. Pate <sup>a,b</sup>






© 2014 BrightSpec, Inc. <sup>a</sup>) 770 Harris St. Suite 104B, Charlottesville, VA 22903, <sup>b</sup>) Department of Chemistry University of Virginia – Charlottesville, VA

### FT-MRR Spectroscopy for Chemical Analysis of SMALL, POLAR, VOLATILES

$$I = MR^2$$

Rotational spectroscopy is high resolution and structure specific. No chemical separation is necessary. The rotational spectrum can be calculated with a high degree of accuracy based on the molecule's mass distribution.

SMALL (less than 150 amu)



There are two effects as the molecule gets heavier:

- 1) The rotational energy levels get closer, spreading the population thinner (weaker emission)
- 2) Conformers and vibrationally excited modes contribute distinct spectra (weaker emission, dense spectra)


### POLAR



Stronger coupling between the electric field of the excitation light and the oscillating dipole means stronger signals.

### VOLATILE

High resolution MRR spectroscopy require freely rotating molecules at reduced pressure ( $10^{-4}$  atm). Only 1mL of STP gas is required.



**Liquids:**  
Headspace Solutions  
Reaction monitoring



**Solids:**  
Residual Solvent Identification  
TGA



**Gases:**  
Cylinder purity  
Process control, feed gas  
Airborne contaminants


### FT-MRR New Spectroscopy for Analytical Applications FUNDAMENTAL ADVANTAGES

<b>Selectivity</b>	MRR is high resolution, no complex chemometrics required	<b>Direct Analysis</b> No chromatography
<b>Structure Specificity</b>	MRR spectra are unique to the moments of inertia of a rotating molecule. Any mass re-distribution shifts the spectrum.	<b>Isomers, Conformers, Isotopologues, Enantiomers</b>
<b>Sensitivity</b>	The chirped-pulse Fourier transform technique maximizes the measurement duty cycle and excitation efficiency.	<b>In Minutes: f mols, 100s ppb</b> <b>10<sup>5</sup> Dynamic Range</b>
<b>Precision</b>	Very stable millimeter wave light sources produce consecutive spectra that are reproducible to < 0.5%	<b>On line monitoring</b> Method transfer
<b>Simplicity</b>	No lasers or complicated optics. The highly stable light field is digitally generated, detected, and analyzed	<b>Direct Automated</b> Easy Operation
<b>Multi-Species Sensing</b>	Excitation and detection are digitally controlled for full-band and targeted monitoring of multiple species in a mixture	<b>One instrument, many molecules</b> Reprogrammable
<b>Breadth of Application</b>	MRR is a non-destructive, low pressure gas measurement ( $10^{-4}$ atm)	<b>Solvents</b> Airborne Toxins Gas Mixtures

### Spectroscopy AND Sensitivity

Molecules of Interest	Predicted Detection Limits		
	Detection Limit (40 seconds, pmol)	Residual Solvent Class 3	Detection Limit (40 seconds, pmol)
Acetonitrile*	0.02	Formic Acid	0.28
Carbonyl Sulfide*	0.04	Dimethyl sulfoxide	0.41
Ethylene oxide*	0.05	Acetone*	1.5
Trifluoromethane	0.09	Ethanol*	1.6
Acrylonitrile*	0.11	Acetic Acid	2.8
Propionitrile*	0.14	Methylethyl ketone	3.7
Acetaldehyde*	0.20	2-propanol*	4.2
Nitric Acid	0.16	Ethyl formate	4.5
Formaldehyde (in water)*	0.26	1-propanol*	10
Propyne*	0.33	2-butanol	11
Methyl isocyanate*	0.38	Ethyl acetate	23
Chloroacetonitrile*	0.66	1-butanol	32
Methanesulfonyl chloride*	0.71	Diethyl ether	62
Nitrous Oxide*	1.26	Anisole	80
1-Butyne*	2.9	3-methyl-1-butanol	140
Nitrogen Dioxide	8.6	tert-butylmethyl ether	370
Methyl acetate	11	Methylisobutyl ketone	372

\* Experimentally measured at BrightSpec Labs






### Composition of solution headspace

Experimental:

- 1) A standard sample of each component of VOC Mix 6 dissolved in methanol was diluted 1:10 in water. (20 – 500 ug/mL)
- 2) A standard sample of EPA VOC Mix 6 was diluted to 1:20 in water. (each component 100 ug/mL)
- 3) 1mL of each solution was injected through a septum into an evacuated glass vial.
- 4) After 10 minutes of equilibration at room temperature, the FT-MRR spectrum is recorded for each sample by drawing 10 mTorr of headspace into the spectrometer sample cell via a needle connected to the vacuum transfer line.


Results: Dynamic Range  $10^5$



With 50,000 independent data channels, the high resolution FT-MRR spectrum enables high dynamic range, direct mixture analysis. High sensitivity results can be obtained in 5 minutes for broad searches (above) or in 40 seconds for targeted detection (left). Linearity is retained across a dynamic range of  $10^5$ .

Quantitation in liquid-vapor equilibrium:

1. Reference measurement of A



2. Linearly extrapolate: A in mixture

$$\frac{Signal_{mixture}}{Signal_{reference}} = \text{Partial Pressure } (P_A)$$

3. Apply Henry's Law  $P_A = k_H X_A$

$$X_A = \text{mole fraction A in solution}$$

Solvent Detection Limits in Water

	Broadband (5min)		Targeted (40sec)
Chloromethane	0.005 ug/mL	0.1 ug/L	0.1 ug/L
Bromomethane	0.022 ug/mL	0.8 ug/L	0.8 ug/L
Chloroethane	0.026 ug/mL	0.5 ug/L	0.5 ug/L
Vinyl Chloride	0.042 ug/mL	0.5 ug/L	0.5 ug/L
Dichlorodifluoromethane	100 ug/mL	100 ug/L	100 ug/L
Trifluorochloromethane	100 ug/mL	1000 ug/L	1000 ug/L

For comparison, the EPA method 624 for VOC detection in water lists GC/MS purge and trap detection limits on the order of 0.5 – 1 ug/L. The EPA 624 method requires approximately 60 minutes. The FT-MRR detection technique offers time savings and simplicity. The volatiles here were measured by a direct headspace measurement without heating, trapping, or salt enhancements and with a total measurement time 10 times faster than the established EPA method.


### Implications and Applications

MIXTURE ANALYSIS WITHOUT A CHEMICAL SEPARATION:


- 10 times higher sample throughput
- Reduction in consumables cost (no columns or carrier gases)
- Reduction in resource intensive calibration
- Turn key, automated operation (minimal training)

Vapor phase monitoring is important to many processes in chemical manufacturing. FT-MRR promises new analytical capability enabled by molecular structure selectivity, high sensitivity, high resolution and operational simplicity.

### Formulation: Online Monitoring



### Solutions: Residual Solvent



### Solids: Residual Solvent

