

# No Standards, No Problem: Molecular Rotational Resonance Spectroscopy Cracks the Chiral Code


## Overview

The proportion of (R)- and (S)-enantiomers in an essential oil directly influences its aroma, biological activity, and therapeutic properties. For researchers in this industry, accurately determining the enantiomeric composition of a sample is critical, as the ratio of (R)- to (S)-enantiomers provides a unique fingerprint for authenticity, quality control, and batch consistency. Traditional methods for chiral analysis, such as chiral gas chromatography (GC) and high-performance liquid chromatography (HPLC), rely on chiral stationary phases or derivatization agents to separate enantiomers. Although effective, these techniques require extensive sample preparation, lengthy analysis times, and costly enantiopure standards that can be difficult to obtain and verify. Even advanced techniques like chiral GC-MS—despite improved detection capabilities—still rely heavily on these reference materials, adding complexity, cost, and supply chain limitations to the analytical process.

Molecular rotational resonance (MRR) spectroscopy is a compelling alternative for enantioselective analysis, eliminating the need for reference standards while delivering high-resolution chiral differentiation. The first commercially available MRR spectroscopy instrument, the BrightSpec isoMRR™ platform, provides analytical chemists in the essential oils industry with a precise, versatile, and efficient solution for determining (R)/(S) ratios for multiple chiral components. In this application note, we demonstrate the accuracy and reproducibility of this advanced platform in distinguishing between closely related chiral molecules in copaiba, cypress, and nutmeg oil samples.



## Tagging the Truth: MRR for Chiral Authentication of Essential Oils



MRR spectroscopy leverages quantum mechanical predictions to verify the rotational frequencies of chiral complexes, removing the need for costly reference standards. The BrightSpec-MRR platform redefines enantioselective analysis by combining the detailed structural insights you get from nuclear magnetic resonance (NMR) with the speed and throughput needed to support modern analytical chemistry labs. A key advantage of MRR is its chiral tagging capability, which allows for the precise differentiation of enantiomers. This technique introduces a small chiral molecule, such as propylene oxide (PO) (**Figure 1**), to form non-covalent diastereomeric complexes with analytes. These diastereomeric complexes produce distinct rotational spectra, enabling accurate determination of absolute configuration (R or S) and enantiomeric excess (EE%) without relying on enantiopure reference standards. By directly measuring enantiomeric composition, MRR simplifies chiral analysis, offering analytical chemists a more efficient and accessible method for product authentication and quality control.

## Using MRR Spectroscopy to Distinguish Chiral Signatures of Essential Oils

### Copaiba Oil

$\beta$ -caryophyllene is notable among terpenes because it can interact with the body's endocannabinoid system (CB2 receptors), contributing to its reported anti-inflammatory and therapeutic effects.<sup>1-3</sup> As a major component of copaiba oil, it has garnered much interest in the pharmaceutical, cosmetic, and aromatherapy industries. However, its enantiomeric composition can vary naturally, making accurate enantioselective analysis essential for product authentication and quality control in commercial applications.

The data presented in **Table 1** demonstrate the effectiveness of the BrightSpec-MRR platform for efficiently and precisely determining the enantiomeric composition of  $\beta$ -caryophyllene-rich essential oils like copaiba. Although  $\beta$ -caryophyllene standards are readily available, enantiopure standards do not exist. **Table 1A** compares the enantiomeric excess values of the  $\beta$ -caryophyllene standard with those measured in copaiba oil to determine whether observed differences might be due to the sample matrix rather than instrumental variability.

Essential Grade Copaiba Oil			
Target Molecule	Major Enantiomer	EE (%)	
$\beta$ -Caryophyllene Standard	S	>96.01	
$\beta$ -Caryophyllene Standard	S	>96.08	
$\beta$ -Caryophyllene Standard	S	>95.69	
Copaiba Oil (Vendor 1)	S	>93.58	
Copaiba Oil (Vendor 1)	S	>93.34	
Copaiba Oil (Vendor 1)	S	>91.74	

High $\beta$ -Caryophyllene Essential Oils (% EE Method)	
Step	Major Enantiomer (S)
MRR Measurement	40 min*
Full Analysis Cycle Time	59 min

\* Depends on the required analysis sensitivity

**Table 1. EE analysis of  $\beta$ -caryophyllene-rich essential grade copaiba oil using MRR. A.** Enantioselective analysis of essential grade copaiba oil and  $\beta$ -caryophyllene standard. The (S)-enantiomer was consistently dominant across all samples. Samples were run in triplicate and obtained from the same vendor. **B.** MRR measurement and full analysis cycle time for determining EE in high  $\beta$ -caryophyllene essential oils. The major enantiomer detected was the (S)-enantiomer, with an MRR measurement time of 40 minutes and a full analysis cycle time of 59 minutes.

The  $\beta$ -caryophyllene standard consistently exhibited a strong (S)-enantiomer preference, with EE values greater than 96% (**Table 1A**). Similarly, essential-grade copaiba oil samples also showed a dominant (S)-enantiomer preference, with EE values greater than 92% across triplicate analyses. MRR spectroscopy offers a practical alternative to traditional chiral separation methods. With an MRR measurement completed in just 40 minutes and a full analysis cycle under 60 minutes (**Table 1B**), this approach considerably reduces the analysis time and complexity associated with traditional methods like chiral GC or HPLC, which typically require longer run times and extensive sample preparation.

### Cypress Oil

Cypress oil, particularly its major component  $\alpha$ -pinene, is a popular ingredient in respiratory blends for its bronchodilator effects and congestion relief. The (R)-enantiomers of key monoterpenes like  $\alpha$ -pinene and 3-carene, which dominate cypress oil, contribute to its biological activity and consumer appeal. Using the BrightSpec-MRR platform, we performed enantioselective analysis on essential grade cypress oil samples, as shown in **Table 2**, revealing distinct enantiomeric distributions for these two terpenes.

(R)- $\alpha$ -pinene emerged as the dominant enantiomer across all samples; however, its EE varied moderately, ranging from 22.76% to 35.08% in triplicate measurements. This lower enantiopurity indicates that  $\alpha$ -pinene's chiral composition in cypress oil is less consistently enantioenriched, potentially reflecting influences from biosynthetic pathways, environmental conditions, or sample matrix interactions. In contrast, (R)-3-carene displayed a consistently high enantiomeric excess (EE > 91%), demonstrating stronger and more consistent enantioselectivity during biosynthesis (**Table 2A**).

A.

Essential Grade Cypress Oil		
Target Molecule (Vendor 2)	Major Enantiomer	EE (%)
$\alpha$ -pinene	R	22.76
$\alpha$ -pinene	R	35.08
$\alpha$ -pinene	R	29.69
3-carene	R	>92.29
3-carene	R	>93.03
3-carene	R	>91.44

B.

Two Component Essential Oils (% EE Method)	
Step	Major Enantiomer (R)
MRR Measurement	25 min*
Full Analysis Cycle Time	44 min

\* Depends on the required analysis sensitivity

**Table 2. EE analysis of essential grade cypress oil using MRR. A.** Enantioselective analysis of essential grade cypress oil monoterpenes  $\alpha$ -pinene and 3-carene. The (R)-enantiomer was consistently dominant across all samples. Samples were run in triplicate and obtained from the same vendor. **B.** MRR measurement and full analysis cycle time for determining EE in two-component essential oils. The major enantiomer detected was the (R)-enantiomer, with an MRR measurement time of 25 minutes and a full analysis cycle time of 44 minutes.

Together, these findings illustrate that although both  $\alpha$ -pinene and 3-carene predominantly exist as (R)-enantiomers, 3-carene is considerably more enantiopure. The distinct chiral profiles identified through MRR spectroscopy thus provide a valuable chemical fingerprint for verifying the authenticity and quality of cypress oil, highlighting the platform's efficiency (**Table 2B**) and precision in rapid, enantioselective analysis.

## Nutmeg Oil

Nutmeg oil is well known for its calming, warming, and pain-relieving properties, often utilized to reduce stress, alleviate muscular discomfort, and promote digestive wellness. Monoterpenes such as  $\alpha$ -pinene and  $\beta$ -pinene, commonly found in nutmeg and other essential oils, have distinct structural characteristics, aromas, and therapeutic properties. Clearly understanding these differences helps formulators select essential oils with the ideal profiles tailored to specific therapeutic or aromatic needs. In our analysis of a provided nutmeg oil sample, (S)- $\alpha$ -pinene was the predominant enantiomer identified, present with moderate EE (approximately 47–51%), indicating a balanced chiral composition that is characteristic of authentic nutmeg essential oils. This balanced chiral composition reflects the natural variability typical of authentic nutmeg essential oils, highlighting the importance of precise analytical methods like MRR spectroscopy for verifying product authenticity and quality control.

In contrast, the much lower EE values observed for  $\beta$ -pinene (approximately 0–4%) suggest an almost racemic mixture, pointing to minimal enantiomeric selectivity during its biosynthesis in nutmeg. From an aromatherapeutic standpoint, the balanced enantiomeric profile of nutmeg oil—particularly the notable presence of (S)- $\alpha$ -pinene—contributes to its

A.

Essential Grade Nutmeg Oil			
Target Molecule (Vendor 2)	Major Enantiomer	EE (%)	Uncertainty
$\alpha$ -pinene	S	50.59	7.28
$\alpha$ -pinene	S	47.38	7.58
$\alpha$ -pinene	S	50.76	7.27
$\beta$ -pinene	R	0.16	4.65
$\beta$ -pinene	S	3.98	4.65
$\beta$ -pinene	S	3.19	4.65

B.

Three Component Essential Oil (% EE Method)	
Step	Major Enantiomer (R)
MRR Measurement	25 min*
Full Analysis Cycle Time	44 min

\* Depends on the required analysis sensitivity

**Table 3. EE analysis of essential grade nutmeg oil using MRR. A.** Enantioselective analysis of essential grade nutmeg oil monoterpenes (S)- $\alpha$ -pinene and  $\beta$ -pinene. (S)- $\alpha$ -pinene was the predominant enantiomer. Samples were run in triplicate and obtained from the same vendor. **B.** MRR measurement and full analysis cycle time for determining EE in three-component essential oils. The major enantiomer detected was the (R)-enantiomer, with an MRR measurement time of 25 minutes and a full analysis cycle time of 44 minutes.

characteristic spicy, warming aroma. The moderate uncertainty values reported (3.91–7.58) further confirm MRR spectroscopy as a robust analytical technique capable of accurately quantifying enantiomeric distributions, reinforcing consumer confidence in nutmeg oil's therapeutic and aromatherapeutic applications.

## Conclusion

MRR spectroscopy, as implemented on BrightSpec's novel isoMRR platform, represents a transformative advancement in the chiral analysis of essential oils and fragrances. By eliminating the need for enantiopure standards, MRR substantially reduces costs, simplifies analytical processes, and mitigates challenges associated with sourcing reference materials. MRR can also measure multiple chiral components out of the same oil, avoiding the costly process of screening columns, and the extensive method development needed in chiral chromatography.

Its capability to deliver rapid and high-resolution differentiation of closely related chiral molecules enables essential oil producers to confidently authenticate products, ensure batch-to-batch consistency, and meet quality guidelines. The detailed analyses of copaiba, cypress, and nutmeg oils highlight the practical utility of MRR for accurately determining enantiomeric compositions across complex matrices. By leveraging chiral tagging and precise spectral measurements, MRR provides unmatched specificity, reliability, and efficiency compared to traditional methods such as GC and HPLC. Ultimately, this innovative approach reinforces consumer confidence and satisfaction by consistently validating the authenticity of essential oil products.

## References

1. Klauke, A.-L., Racz, I., Pradier, B., et al. The Cannabinoid CB<sub>2</sub> Receptor-Selective Phytocannabinoid  $\beta$ -Caryophyllene Exerts Analgesic Effects in Mouse Models of Inflammatory and Neuropathic Pain. *Eur. Neuropsychopharmacol.* 2014, 24;608–620.
2. Gertsch, J., Leonti, M., Raduner, S., et al. Beta-Caryophyllene Is a Dietary Cannabinoid. *Proc. Natl. Acad. Sci. U.S.A.* 2008, 105;9099–9104.
3. Fidy, K., Fiedorowicz, A., Strz̄dała, L., et al.  $\beta$ -Caryophyllene and  $\beta$ -Caryophyllene Oxide—Natural Compounds of Anticancer and Analgesic Properties. *Cancer Med.* 2016, 5;3007–3017.
4. Sonstrom, R. E., Cannon, D. M., Neill, J. L. Chiral Analysis of Linalool, an Important Natural Fragrance and Flavor Compound, by Molecular Rotational Resonance Spectroscopy. *Symmetry* 2022, 14 (509).
5. Wilson, T. M., Davis, A., Sonstrom, R. E., et al. Essential Oil Composition and Enantioselective Profile of *Agastache urticifolia* (Lamiaceae) and *Monardella odoratissima* (Lamiaceae) from Utah. *Molecules* 2023, 28 (2249).
6. Wang, M., Lee, J., Zhao, J., et al. Comprehensive Quality Assessment of Peppermint Oils and Commercial Products: An Integrated Approach Involving Conventional and Chiral GC/MS Coupled with Chemometrics. *J. Chromatogr. B* 2024, 1232, 123953. <https://doi.org/10.1016/j.jchromb.2023.123953>.