

High-Throughput, Fully Automated Analysis of USP<467> Residual Solvents Using Headspace-MRR

Alexander V. Mikhonin, Dreka Burgess, Voislav Blagojevic, and Justin L. Neill
BrightSpec, Inc., Charlottesville, VA

Introduction

USP <467> recommends static headspace GC (SH-GC) for residual solvent analysis. Yet, many important solvents, including low-volatility compounds, acids and amines remain analytically challenging.

Molecular Rotational Resonance (MRR) spectroscopy directly identifies volatile molecules in mixtures by their unique gas-phase rotational spectra, providing unambiguous ID and reliable quantitation without chromatography.

Here we demonstrate fully automated headspace-MRR method for high-throughput analysis of GC-unfriendly residual solvents in pharmaceutical substances, expanding analytical coverage beyond conventional SH-GC.

MRR Applications

Structure verification when standards are not available

- J.L. Neill et al., *Org. Process Res. Dev.*, **2019**, 23 (5), 1046.
- J.L. Neill et al., *Anal Sci Adv.* **2023**, 4, 204

Direct analysis of stereo- and regio-isomer mixtures

- J.L. Neill et al., *J. Pharm. Biomed. Analysis*, **2020**, 189, 113474
- L.A. Joyce et al., *Chem. Sci.*, **2020**, 11, 6332
- R.E. Sonstrom et al., *Chirality*, **2022**, 34 (1), 114
- J.L. Neill et al., *Anal Sci Adv.* **2023**, 4, 204

Online, highly-selective reaction monitoring

- J.L. Neill et al., *Org. Process Res. Dev.*, **2019**, 23 (5), 1046.
- A.A. Byars et al., *Precis. Chem.*, **2024**, 2, 57

Site-specific isotope analysis, deuterated isotopomers

- Z.P. Vang et al., *J. Am. Chem. Soc.*, **2021**, 143 (20), 7707
- M.D. Mills et al., *Angew. Chemie*, **2022**, 61 (33), e202207275
- R.E. Sonstrom et al., *Org. Process Res. Dev.*, **2023**, 27 (7), 1185

GC-Unfriendly Volatile Impurities


- A.V. Mikhonin et al., *USP Stimuli Article*, **USP-NF**, March 3, 2025
- Formaldehyde, ethylene oxide and acetaldehyde in pharmaceutical substances
- Ethylene glycol and diethylene glycol in high-risk excipients

Solvent Composition / Solvent Exchange

- A.V. Mikhonin et al., "Direct Solvent Composition Monitoring Using Headspace-MRR: Implications for Crystallization Process Optimization", Poster at PittCon 2025

Today's focus


Headspace-MRR Workflow



Samples per batch: up to 66
Time to first results: about 1 hour
Analysis cycle time: 3 to 30 minutes per sample


Headspace-MRR Features

- Spectroscopy, direct analysis without separation
- Unique molecular fingerprints, no spectral overlaps
- Reliable, highly-selective quantitation in mixtures
- Sensitive: no limits on headspace volumes




- Bridges analytical gaps of conventional methods
- Improves sample throughput

MRR Selectivity



MRR Method Validation



Why MRR for Impurities



GC-Unfriendly Impurities

- ✓ High-boilers from USP<467> Class 2 Mixture C solvents
- ✓ EG and DEG in excipients
- ✓ Small, highly-volatile species: formaldehyde, ethylene oxide, dimethyl sulfide, etc.
- ✓ Water-soluble acids
- ✓ Nitrosamines

Fast Analysis, Simple Workflow

- ✓ No derivatization
- ✓ No wait for separation
- ✓ 1 to 10 minutes per sample
- ✓ Analytes from different standard mixtures in one run

Low Operational Costs

- ✓ No mobile-phase solvents
- ✓ Fully automated "push-button" analysis

Structure ID without Standards

- ✓ Quasi-unknowns: unambiguous "yes" or "no" to a candidate structure (volatile NDSRI, for example)

References

A.V. Mikhonin et al., "Direct Analysis of Class 2 Residual Solvents Using Molecular Rotational Resonance Spectroscopy", *USP Pharmacopeial Forum*, PF 51(2), March 3, 2025.

MRR Analysis Metrics

USP<467> Class 2 Mixture C Residual Solvent	2-Methoxyethanol	2-Ethoxyethanol	Dimethylformamide	Dimethylacetamide	N-Methylpyrrolidone	
USP<467> Limit in substance (μg/g)	50	160	880	1090	530	USP Pass Criteria
Limit in 20 mg/mL substance solution (μg/mL)	1	3.2	17.6	21.8	10.6	
Confirmed MRR Linear Range in Solution (μg/mL)	0.33 – 2.5	0.51 – 8	3.03 – 44	2.27 – 54.5	3.39 – 26.5	NLT 50-150% of USP limit
Linearity (R ²)	0.97	0.98	0.98	0.99	0.99	R ² NLT 0.9
Repeatability (%)	8.7%	9.8%	6.2%	9.1%	11.6%	NMT 20%
Recovery in 20 mg/mL Acetaminophen (%)	96.8%	92.3%	100.8%	102.7%	109.7%	80% – 120%
Low Limit of Quantitation, LOQ, in Solution (μg/mL)	0.33	0.51	3.03	2.27	3.39	NMT 50% of USP Limit
LOQ in 20 mg/mL Acetaminophen (μg/g)	16	26	152	113	169	
MRR Metrics vs USP<467> and USP<1467> Requirements	pass	pass	pass	pass	pass	pass

Conclusions

- Reliable ID and quantitation with one instrument
- MRR can bridge analytical gaps of conventional methods
 - ❖ MRR can help with selectivity
 - ❖ MRR can help with sensitivity (use-case: low volatile analytes in injection-unfriendly matrices)
- MRR can help with analysis throughput